



Data Sheet (DS) 7TM0040N

Issuing Date: 2026-07-29

Description	BB3 (non-phospho), Bombesin Receptor 3 Antibody
Format	Purified, Liquid
Product Type	Rabbit Polyclonal Antibody
Isotype	Polyclonal IgG
Quantity	100 µl

Product Details

Applications	This product has been reported to work in the following applications:	
		Dilution
	Western Blot	1:1000
	Immunohistochemistry	1:100
	Immunocytochemistry	-
	This information is derived from testing within our laboratories and peer-reviewed publications. Please refer to references indicated for further information. For general protocol recommendations, please visit https://7tmantibodies.com/7tm-antibodies-support/7tm-protocols/ Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.	
Target Species	Human	
Product Form	Purified IgG, liquid	
Antiserum Preparation	Antiserum to BB3 was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.	
Immunogen	A synthetic peptide with sequence EISVTSFTGCSVKQAEDRF corresponds to amino acids 381-399 of human BB3.	
Storage Buffer	Dulbecco's PBS, pH 7.4, with 150 mM NaCl, 0.02% sodium azide	
Specificity	The BB3 receptor antibody is directed against the distal end of the carboxyl-terminal tail of human Bombesin Receptor 3. It can be used to detect BB3 in formalin-fixed, paraffin-embedded tissue sections by immunohistochemistry. It can be used to detect total BB3 receptors in Western blots independent of phosphorylation. The BB3 antibody can also be used to isolate and enrich BB3 receptors from cell and tissue lysates.	
Guarantee	12 months from date of dispatch	
Storage	Store at -20°C.	

	This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.
Regulatory	For research purposes only
Health and Safety Information	Material Safety Data Sheet documentation is available at https://7tmantibodies.com/phosphosite-7tm-antibodies/bombesin-receptors/bb3brs3/802/bb3-non-phospho-bombesin-receptor-3-antibody in the downloads section as: Safety Data Sheet EU Safety Data Sheet US

Related Products

[pT360/pS361-BB3 \(Bombesin Receptor 2/Gastrin-Releasing Peptide Receptor Antibody\)](#)

[pT363/pT364-BB3 \(Bombesin Receptor 2/Gastrin-Releasing Peptide Receptor Antibody\)](#)

[pT370-BB3 \(Bombesin Receptor 2/Gastrin-Releasing Peptide Receptor Antibody\)](#)

Details of the Supplier of the Data Sheet

Supplier

7TM Antibodies GmbH

Hans-Knöll-Str. 6

07745 Jena – Germany

TEL: ++49 151 20130575

FAX: ++49 3641 241 49 58

Email: info@7tmantibodies.com

Website: 7tmantibodies.com